



Recombinant Nuclear hormone receptor family member nhr-136 (nhr-136)

Product Code	CSB-YP518838CXY
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	O01930
Product Type	Recombinant Protein
Immunogen Species	Caenorhabditis elegans
Purity	>85% (SDS-PAGE)
Sequence	MLPDFITVPS TSEKCLSP EL LLPNMDENLE DSKPSSLKYQ KAMMSRGTCP SNCKVCRHSA TGYHYDVPSC NGCKTFFRRS ILDGRKYTCL KMRKCLSGTE PVDLSRRMCR ACRFEKCVEA GMNPSAIQAD MKTTD GELLR KEIMIKQKTA VDFLNT PQVI MSFEDKVQGI IGKLTVMELK IEPLYTGGLP PGNRDIRKLD ELIDAPLILS YDEIPNLKYC PSVDELTGEI KPSGACYIHC GYLASIEYSK MFDFAHKIDV ASKATLIKHA TIMCADIMTA FFSYYQRKSD RLIHPNGMFA GPPKYRYGEA GTKYQASMQR TLATVLRHEL NRIEYMLLKA IVLCNPAVSS LSISVQQIIG KEREYVRTL LTYCLLN YGS VHGPSRFSAL LAIMSVLESQ QKNAKDFHLL AKATILKDAV RYTRISNLYE QIMES
Source	Yeast
Target Names	nhr-136
Protein Names	Recommended name: Nuclear hormone receptor family member nhr-136
Expression Region	1-435
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.