



Recombinant Arabidopsis thaliana Glutathione S-transferase U21 (GSTU21)

Product Code	CSB-MP518988DOA
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	F4IA73
Product Type	Recombinant Protein
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Purity	>85% (SDS-PAGE)
Sequence	MAAEVILLGF WPSMFGMRTM IALEEKGVKY EYREEDVINN KSPLLLEMNP IHKTIPVLIH NGKPVLESLI QIQYIDEVWS DNNSFLPSDP YHRAQALFWA DFIDKKEQLY VCGRKTWATK GEELEAANKE FIEILKTLQC ELGEKPYFGG DKFGFVDIVL IGFYSWFPAY QKFGNFSIEP ECLKLIAWVK RCMQRESVAK ALPDSEKVVG YVLQLKKLYG IE
Source	Mammalian cell
Target Names	GSTU21
Protein Names	Recommended name: Glutathione S-transferase U21 Short name= AtGSTU21 EC= 2.5.1.18 Alternative name(s): GST class-tau member 21
Expression Region	1-222
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.