



Recombinant Pneumococcus phage Dp-1 Lysin (PAL)

Product Code	CSB-MP519054EXF
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	O03979
Product Type	Recombinant Protein
Immunogen Species	Pneumococcus phage Dp-1 (Bacteriophage Dp-1)
Purity	>85% (SDS-PAGE)
Sequence	MGVDIEKGVA WMQARKGRVS YSMDFRDGPD SYDCSSSMYY ALRSAGASSA GWAVNTEYMH AWLIENGYEL ISENAPWDAK RGDIFIWGRK GASAGAGGHT GMFIDSDNII HCNAYYDGIS VNDHDERWYY AGQPYYYVYR LTNANAQPAE KKLGWQKDAT GFWYARANGT YPKDEFEYIE ENKSWFYFDD QGYMLAEKWL KHTDGNWYWF DRDGYMATSW KRIGESWYYF NRDGSMTGW IKYYDNWYYC DATNGDMKSN AFIRYNDGWY LLLPDGRLAD KPQFTVEPDG LITAKV
Source	Mammalian cell
Target Names	PAL
Protein Names	Recommended name: Lysin EC= 3.5.1.28 Alternative name(s): Cell wall hydrolase Lytic amidase N-acetylmuramoyl-L-alanine amidase
Expression Region	1-296
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.