



Recombinant Schizosaccharomyces pombe Protein disulfide-isomerase C17H9.14c (SPAC17H9.14c)

Product Code	CSB-EP519308SXV
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	O13811
Product Type	Recombinant Protein
Immunogen Species	Schizosaccharomyces pombe (strain 972 / ATCC 24843) (Fission yeast)
Purity	≥85% (SDS-PAGE)
Sequence	S GVVELQSLNE LENTIRASKK GALIEFYATW CGHCKSLAPV YEELGALFED HNDVLIGKID ADTHSDVADK YHITGFPTLI WFPPDGSEPV QYSNARDVDS LTQFVSEKTG IKKRKIVLPS NVVELDSLNF DKVVMDDKKD VLVEFYADWC GYCKRLAPTY ETLGKVFKNE PNVEIVKINA DVFADIGRLH EVASFPTIKF FPKDDKDKPE LYEGDRSLES LIEYINKKSG TQRSPDG TLL STAGRIPTFD EFAAEFLDMS NAAKEVVLEK VKQLALEDSS RWTKYYYKKVF EKILNDENWV HKEAKRLSKL LRQKSIALAS ADDFKTRLNI LNSFLPGNH
Source	E.coli
Target Names	SPAC17H9.14c
Protein Names	Recommended name: Protein disulfide-isomerase C17H9.14c EC= 5.3.4.1
Expression Region	20-359
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.