



Recombinant Schizosaccharomyces pombe Uncharacterized protein C1A6.10 (SPAC1A6.10, SPAC30D11.05c)

Product Code	CSB-MP519315SXV
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	O13861
Product Type	Recombinant Protein
Immunogen Species	Schizosaccharomyces pombe (strain 972 / ATCC 24843) (Fission yeast)
Purity	>85% (SDS-PAGE)
Sequence	METLQVEVFF LFDFTSIMAG FKVPSWITYK SFWIAVSSSV TTACVILGTL EFRKHRSIRR LQSMIVPEAG KSIQLSSSGV PIEIYDAGEE DEGISKGV DENLIREQLA RNYAFFGEDG MERLRNSFVI VVGCGGVGSW VINMLARSGV QKIRIVDFDQ VSLSSLNRHS IATLQDVGTP KTLAIKKAIA KFAPWIEVDA RNALFNPDSA DDLLSGNPFD VIDAIDNIQT KVDLLSYCYN HKLPVIASTG SACKSDPTRV NIADISATSE DPLSRATRRR LRLLGIMEGI PVVFSTKPD PRKASLLPLS EEEFEKGDVD ELSALPEFRA RILPVIGPMP GIFGLTIATY VLTSIAKYPM DPISTMTRPR LYEEAVKRLH AEARKAGVNL DKTFNASEMS YIIIEVYVGR SALPPHESQK VTVVRWNPQL PFDHTNLVAM TRNEARYHED NVLAKNVDPS TVYGKDVIEV VHSFLRRLRM WEMLY
Source	Mammalian cell
Target Names	SPAC1A6.10
Protein Names	Recommended name: Uncharacterized protein C1A6.10
Expression Region	1-485
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.