



Recombinant Schizosaccharomyces pombe Uncharacterized methyltransferase C1B3.06c (SPAC1B3.06c)

Product Code	CSB-EP519317SXV
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	O13871
Product Type	Recombinant Protein
Immunogen Species	Schizosaccharomyces pombe (strain 972 / ATCC 24843) (Fission yeast)
Purity	>85% (SDS-PAGE)
Sequence	MSKPTDQLYY ANGVDSYIAE THAWRTPETC STYMLKYVKK TDRILDVGC PGTITVGFPK YVPEGEVIGV EPSQELLDKA EEALRKEETL KKEKINNC RLGSIYKLPF PDNTFDIVNT HQVLVHLQDP VAALVELKRV TKPGGYV EADLLSACVY PKEYEHDLLL QSQARINLHG TNPTAGRSLR GWAI DAKYVA ENIHSSASTW CFADEETRKW VSRLFIQRLV HSNERLDDDD AKDQSLRKRVAEAWQRWKED SRGCFMFDG QIVYKKEE
Source	E.coli
Target Names	SPAC1B3.06c
Protein Names	Recommended name: Uncharacterized methyltransferase C1B3.06c EC=2.1.1.-
Expression Region	1-278
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.