



# Recombinant *Schizosaccharomyces pombe* Probable ATP-citrate synthase subunit 2 (SPAC22A12.16)

<b>Product Code</b>	CSB-BP519321SXV
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	O13907
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	<i>Schizosaccharomyces pombe</i> (strain 972 / ATCC 24843) (Fission yeast)
<b>Purity</b>	≥85% (SDS-PAGE)
<b>Sequence</b>	MSAKSIREYD GKAVLAYWLN RSPSISKEEY KTVSATPAVQ LAQIQFPLPT LVIPAESTTY REVVEDVFAK VEQEHPWVKE TKLVAKPDQL IKRRGKSGLL KLNATWDEAK EWIRERAGKN QKVQHAVGYL TFLVPEFVP HPPNTEYYIN INSVREGDWI LFCNEGGVDV GDVDAKARKL LVPVRLSEFP SRATIASTLL SDIPVEQHEs LVDFIIRLYS VYVDCQFTYL EINPLVVIPT AKGADV FYLD LAAKLDQTAE FECGAKWAVA RAPESLGIKT SGEESGAINA DHGPPMV FPA PFGRELSKEE AYVQGLDAKT GASLKLTLN AEGRVWNLVA GGGASVVYAD AVAVNGAADE LANYGEYS GA PTDGQTYEYA KTVLDLMTRG EPRADGKVL F IGGGIANFTS PAVTFRAIAR ALGDYKDKLH AHKVSIVVRR AGPNYQEGLR VIREAGKKFD LPLKVY GPEC HISGIVPMGL GKAPVEAEWQ MA
<b>Source</b>	Baculovirus
<b>Target Names</b>	SPAC22A12.16
<b>Protein Names</b>	Recommended name: Probable ATP-citrate synthase subunit 2 EC= 2.3.3.8 Alternative name(s): ATP-citrate (pro-S-)-lyase 2 Citrate cleavage enzyme subunit 2
<b>Expression Region</b>	1-492
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients,



storage temperature and the stability of the protein itself.  
Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.