



Recombinant *Schizosaccharomyces pombe* Probable ATP-citrate synthase subunit 2 (SPAC22A12.16)

Product Code	CSB-EP519321SXV-B
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	O13907
Product Type	Recombinant Protein
Immunogen Species	<i>Schizosaccharomyces pombe</i> (strain 972 / ATCC 24843) (Fission yeast)
Purity	≥85% (SDS-PAGE)
Sequence	MSAKSIREYD GKAVLAYWLN RSPSISKEEY KTVSATPAVQ LAQIQFPLPT LVIPAESTTY REVVEDVFAK VEQEHPWVKE TKLVAKPDQL IKRRGKSGLL KLNATWDEAK EWIRERAGKN QKVQHAVGYL TFFLVEPFVP HPPNTEYYIN INSVREGDWI LFCNEGGVDV GDVDAKARKL LVPVRLSEFP SRATIASTLL SDIPVEQHEs LVDFIIRLYS VYVDCQFTYL EINPLVVIPT AKGADVIFYLD LAAKLDQTAE FECGAKWAVA RAPESLGIKT SGEESGAINA DHGPPMVFPA PFGRELSKEE AYVQGLDAKT GASLKLTLN AEGRVWNLVA GGGASVVYAD AVAVNGAADE LANYGEYSGA PTDGQTYEYA KTVLDLMTRG EPRADGKVLV IGGGIANFTS PAVTFRAIAR ALGDYKDKLH AHKVSIVVRR AGPNYQEGLR VIREAGKKFD LPLKVYGPPEC HISGIVPMGL GKAPVEAEWQ MA
Source	E.coli
Target Names	SPAC22A12.16
Protein Names	Recommended name: Probable ATP-citrate synthase subunit 2 EC= 2.3.3.8 Alternative name(s): ATP-citrate (pro-S-)-lyase 2 Citrate cleavage enzyme subunit 2
Expression Region	1-492
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients,



storage temperature and the stability of the protein itself.
Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.