



Recombinant Schizosaccharomyces pombe Putative ubiquitin thioesterase otu1 (otu1)

Product Code	CSB-YP519333SXV
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	O13974
Product Type	Recombinant Protein
Immunogen Species	Schizosaccharomyces pombe (strain 972 / ATCC 24843) (Fission yeast)
Purity	>85% (SDS-PAGE)
Sequence	MSSLRLRLKY ENQSAVETVE ANATVGSFLD LVAAKFSLPR NSIALKFGFP PQDIPLVNSD VPLSTLVSSG QQILVLKNAA TSFSTNEPAK PPIPNAATKP TFPPQTEISN PPAVSHQSKN TSQDPPYVST PIGDIALRVM PDDNSCLFRA LSKPLGFSPY ELREIVANQV LSNPDIYSTA ILGKPSIEYA SWIRKETSWG GYIELSILSS HFGVEICSVD VKTGRVDSYN PQPATGQRTY IVYSGIHYDL AALAAVLWDT DVDVVLFDAS DVTITPYVQQ LASLLKNMHY YTDTASFSIR CTICGTGLVG EKDATAHALA TGHTQFGEY
Source	Yeast
Target Names	otu1
Protein Names	Recommended name: Putative ubiquitin thioesterase otu1 EC= 3.4.19.12 Alternative name(s): Meiotically up-regulated gene 141 protein OTU domain-containing protein 1
Expression Region	1-329
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.