



Recombinant Schizosaccharomyces pombe Nuclear speckle splicing regulatory protein 1 homolog (SPAC29A4.06c)

Product Code	CSB-EP519339SXV-B
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	O14009
Product Type	Recombinant Protein
Immunogen Species	Schizosaccharomyces pombe (strain 972 / ATCC 24843) (Fission yeast)
Purity	>85% (SDS-PAGE)
Sequence	MSGTGFRYGL NVMKKKKPNE SSNRITFTED DSSSSEQEHA PIPNSFSSQI TAASDASKDD NYDASIYGFD EFYDSMKSAE REQMELKRLE SKDRRPKYME NLIESAKKRK RDLLARERA LLKNELEGD DATEKFVTSS YKKHREEVQK AIEDRKEEDE KSITDTTTGM KDFYASMLDR QEKIHQAAVE GVQNSKKTGA EIGDAMKGQD TLGSDNLILE AKKKGLKLEL NDNNEIVDQR QILTAGLNIR KSSANNSMRD DKKRNHKSSY KRSLSPSTRY HQDRPDKRHG TYSLEEIDKQ RKEFENRQRL QKEKEFQKSR EAALKIHASR NTTETQVQSA RERYLQRKKK AATNP
Source	E.coli
Target Names	SPAC29A4.06c
Protein Names	Recommended name: Nuclear speckle splicing regulatory protein 1 homolog Alternative name(s): Coiled-coil domain-containing protein 55 homolog
Expression Region	1-355
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.