



Recombinant Schizosaccharomyces pombe UPF0664 stress-induced protein C29B12.11c (SPAC29B12.11c)

Product Code	CSB-MP519343SXV
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	O14032
Product Type	Recombinant Protein
Immunogen Species	Schizosaccharomyces pombe (strain 972 / ATCC 24843) (Fission yeast)
Purity	≥85% (SDS-PAGE)
Sequence	MSINWVTYKE GSEDILLE NECMFDCFDG VAISILCKPP SLKSWTCTKG LLCLTNQRLV YIAKDTDCDF KDFQSPVANL KDTKLNQPPF GANYYS GTVM PVPNGGIPCE AEVKLQFNEG GIFNFVEAWN RLIQRFQEVD SVSRVQHLDLP LPPYHRPSSS QDQPPHYEEA VNKS
Source	Mammalian cell
Target Names	SPAC29B12.11c
Protein Names	Recommended name: UPF0664 stress-induced protein C29B12.11c
Expression Region	1-174
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.