



Recombinant *Schizosaccharomyces pombe* Putative tRNA 2'-phosphotransferase

Product Code	CSB-EP519346SXV-B
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	O14045
Product Type	Recombinant Protein
Immunogen Species	<i>Schizosaccharomyces pombe</i> (strain 972 / ATCC 24843) (Fission yeast)
Purity	≥85% (SDS-PAGE)
Sequence	MYKNMNSHEL IEESNNSGTP ATKSSSKPTK KIRPRNDLVH YSKALSKVLR HTAKANGLQI REDGYIEVDS ILKLPQFRGM GMELLSIVK GNDKKRFTME EVEGVLYIRA NQGHSIKAVQ VPMARIDNAS SIPKVVHGTK KELWPVISKQ GLSRMKRNHI HCATGLYGDP GVISGIRKSC TLYIYIDSAK AMQDGVFEYR SENGVILTEG VNGLLSSKYF SRVETSDGEV LLDAKASPKN NRSDESDQSD PESIDPFCDN LQALSMHELE LLEEKHSNFG YSEGIKGM QVAQSGFDDG FKHGSRLGFQ MGKTIGTLKA KLYIFEENEQ MEILKQELDR LQESAEFHIF VANHKEEILK CIREK
Source	E.coli
Target Names	SPAC2C4.12c
Protein Names	Recommended name: Putative tRNA 2'-phosphotransferase EC= 2.7.1.160
Expression Region	1-365
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.