



# Recombinant Schizosaccharomyces pombe Uncharacterized RING finger protein C2F3.16 (SPAC2F3.16)

<b>Product Code</b>	CSB-YP519353SXV
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	O14099
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Schizosaccharomyces pombe (strain 972 / ATCC 24843) (Fission yeast)
<b>Purity</b>	≥85% (SDS-PAGE)
<b>Sequence</b>	MSLSEYLREF AEIFEEAAQE LTSIFELPPD LLKDITEEDV GTDDSCNKKF LEVQKQKEDE EDEEILKGVD LTQQDSVREK IHEIQSMSQL SEKRKALLMQ KMLMSGYLKY RRTHKESDE NQLSSSDLEK TTYDKEQEIL GCSHYMRNCK VQCFDCHEWY TCRHCHNDAC DHVLERPAVE NMLCMICKSV QPAAQYCKYC KNCMGRYYCN KCKLWDDDPN KSSYHCDDCG ICRIGRGLGD DYFHCKTCGL CLPISVFNTH RCIERSTDCN CPICGEYMFN SRERVIFLSC SHPLHQRCHE EYIRTNYRCP TCYKTIINVN SLFRILDMEI ERQPMPYPYN TWISTIRCND CNSRCDTKYH FLGHKCNSCH SYNTCISSY KPLDHPQVSI PRLATAEDAG MRRLMGHSWD NSDEDFNIFG IHSFI
<b>Source</b>	Yeast
<b>Target Names</b>	SPAC2F3.16
<b>Protein Names</b>	Recommended name: Uncharacterized RING finger protein C2F3.16
<b>Expression Region</b>	1-425
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.