



# Recombinant Schizosaccharomyces pombe Sterol 24-C-methyltransferase (SPBC16E9.05)

|                          |                                                                                                                                                                                                                                                                                                                                                                                                                                                         |
|--------------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Product Code</b>      | CSB-BP519388SXV                                                                                                                                                                                                                                                                                                                                                                                                                                         |
| <b>Storage</b>           | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.<br>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.                                                                                                                                                                      |
| <b>Uniprot No.</b>       | O14321                                                                                                                                                                                                                                                                                                                                                                                                                                                  |
| <b>Product Type</b>      | Recombinant Protein                                                                                                                                                                                                                                                                                                                                                                                                                                     |
| <b>Immunogen Species</b> | Schizosaccharomyces pombe (strain 972 / ATCC 24843) (Fission yeast)                                                                                                                                                                                                                                                                                                                                                                                     |
| <b>Purity</b>            | ≥85% (SDS-PAGE)                                                                                                                                                                                                                                                                                                                                                                                                                                         |
| <b>Sequence</b>          | MSSTALLPPN TDQVLSRR LH GKAAEKKTGL AAIASKDVDE QSRKLQEYFE<br>FWDRNHENES EEDRARRIDG YKSVVNSYYD LATDLYEYGW SQSFHFSRFY<br>KGEAFAQSIA RHEHYLAYRM GIKPGSRVLD VGCGVGGP PAR EITEFTGCNL<br>VGLNNNDYQI SRCNNYAVKR NLDKKQVFVK GDFMHMPFED NTFDYVYAIE<br>ATVHAPSLEG VYGEIFRVLK PGGVFGVYEW VMSDDYDSSI PKHREIAYNI<br>EVGDGIPQMV RKCDAVEAIK KVGFNLLEED DLTDHDNPDL PWYYPLTGDI<br>TKCQNIWDVF TVFRTSRLGK LVTRYSVQFL EKIGVAAKGT SKVGDTLAIA<br>QKGLIEGGET HLFTPMFLMI AKKPETDA |
| <b>Source</b>            | Baculovirus                                                                                                                                                                                                                                                                                                                                                                                                                                             |
| <b>Target Names</b>      | erg6                                                                                                                                                                                                                                                                                                                                                                                                                                                    |
| <b>Protein Names</b>     | Recommended name: Sterol 24-C-methyltransferase EC= 2.1.1.41 Alternative name(s): Delta(24)-sterol C-methyltransferase                                                                                                                                                                                                                                                                                                                                  |
| <b>Expression Region</b> | 1-378                                                                                                                                                                                                                                                                                                                                                                                                                                                   |
| <b>Notes</b>             | Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.                                                                                                                                                                                                                                                                                                                                                     |
| <b>Tag Info</b>          | Tag type will be determined during the manufacturing process.                                                                                                                                                                                                                                                                                                                                                                                           |
| <b>Protein Length</b>    | full length protein                                                                                                                                                                                                                                                                                                                                                                                                                                     |
| <b>Reconstitution</b>    | We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.                                             |
| <b>Shelf Life</b>        | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.<br>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.                                                                                                                                                                      |