



Recombinant Nuclear hormone receptor family member nhr-115 (nhr-115)

Product Code	CSB-MP519422CXY
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	O16662
Product Type	Recombinant Protein
Immunogen Species	Caenorhabditis elegans
Purity	>85% (SDS-PAGE)
Sequence	MKKVLFPCQI CGQNSHGTHF GIVSCRACAA FFRRSVNSKW ARKGCLTNFK DKGSCFCKPC RLRKCV EIGM DASKFQYDRD AISAIKHPKI FPSVSYVGR PEFLMFSDSN VTSQKTFIDV QNLVFEVSR Y LDHGCETPIY AENQLKKLTL GFKLMQFDYQ NVKFFDKIGK AEFIDIEYY FLTVAKWIAH FDEFKLDQS LQIKLLQAIW HVWSKIHKCA STAFYRKNP NAKPTQKILR NVCMDRMHVVH KLDTSWMSDY PTEHVTRFML THVYDFKIV ELLKLDPTD VELTFMFAQL CFEYAGKRFQ GEILKITDHF QQVLSNDLHH YYITDQRRER YFQRLTDLMK VNNLIQRSIW ETRPHRELGR VFEISKLEFS HPEMFDDSGF C
Source	Mammalian cell
Target Names	nhr-115
Protein Names	Recommended name: Nuclear hormone receptor family member nhr-115
Expression Region	1-391
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.