



# Recombinant Ataxin-3 homolog (atx-3)

<b>Product Code</b>	CSB-BP519441CXY
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	O17850
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Caenorhabditis elegans
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MSKDDPINSI FFEHQEAALC AQHALNMLLQ DALYKWQDLR DLAIQMDKME QQILGNANPT PGRSENMNES GYFSIQVLEK ALETFSCLKT NIENPAMVDY KNNPLTARAY ICNLRHWFV LRFKGNQWFE LNSVNRGPKL LSDTYVSMFL HQVSSEGYSI FVVQGVLP RS DADDLISLCP VVPPKVTPKK EQKLEKVMTK FFNTVGKRLG GGSGAPPDSQ EEKDLAIAFA MSMETKDGSE VSRSSAEIDE ENLRKAIELS QAPGPSEPAE IPLLTRSRSS TPPGASEPFS NAEQQRDRQ KFLERFEKKK EERNDEK
<b>Source</b>	Baculovirus
<b>Target Names</b>	atx-3
<b>Protein Names</b>	Recommended name: Ataxin-3 homolog EC= 3.4.19.12 Alternative name(s): Machado-Joseph disease-like protein
<b>Expression Region</b>	1-317
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.