



Recombinant RING finger protein 113 homolog (rnf-113)

Product Code	CSB-EP519442CXY-B
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	O17917
Product Type	Recombinant Protein
Immunogen Species	Caenorhabditis elegans
Purity	>85% (SDS-PAGE)
Sequence	MDLFRKPKKR NAPVVRKKES SSDEDQDSEV KDVIQKRRRT NPMVQSTKQL DASTRRADNS SDDSDSDSDN QDIAVATHSF AASGDAGPSG PRDQGATATL EVDTDYSHDA QAQFERVQQQ LKEGVEKDGK ILYKGSALYG AKEAKDTAKG NAASGYNRVG PVRAPQFLRQ TVRWDFAPDI CKDYKETGFC TFGDSCKFVH DRSDYKHGWE IDEEYEAGKY GAEDDANYEI HEGDDTFPED CFICGNPFVD PIVTKCKHYF CTGCALKSFQ KSSKCPICQQ NTENIMNTAK ELLTYLKRKK QQQKQEAQK EEEKSDSDDE KPHECDDHHH HDHEDEPEEP ENDSNVPEAE EKSDEEQEIM MEDVEGLEGG ENDSSESDDDD AEKD
Source	E.coli
Target Names	rnf-113
Protein Names	Recommended name: RING finger protein 113 homolog
Expression Region	1-384
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.