



Recombinant Nuclear hormone receptor family member nhr-218 (nhr-218)

Product Code	CSB-MP519445CXY
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	O18086
Product Type	Recombinant Protein
Immunogen Species	Caenorhabditis elegans
Purity	>85% (SDS-PAGE)
Sequence	MVASSFQHHL EPMLTQPIPC QICTYQSHGV NFNVMTCRAC AAFFRRSLVC GMRYHCKTRK NDCRIDSTER HFCRLCRFQK CLQMGMKAEK IQQNRDPISS TFPGTSTEPE LSEIVDPENE KSYIFHGTLY GFKSLLAEVR YIFSISRNY DSPLTDLENG FKLITRHQKR RYIDIEDRIN FQNLTDLFR LG HIKNCATWLT HSSFFQSLTE TENLLILKST WHVWSWLELL SVSVEIFGNQ VCEEKIVFLS EKIAVDIVKV FRYILKPLNK QEKRKVEKEL NPIFHILFDD VARKLQNLKP SSLEINYMLW QLVWFVAEKV LNEDNLRHGE QYTNQLASDL HNHYKNDLHL EQYAQRVLKM MAIVKSLQKH LMNIHKIIDY ADNFCNLFAM K
Source	Mammalian cell
Target Names	nhr-218
Protein Names	Recommended name: Nuclear hormone receptor family member nhr-218
Expression Region	1-391
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.