



Recombinant Arabidopsis thaliana Probable glutathione peroxidase 3, mitochondrial (GPX3)

Product Code	CSB-MP519517DOA
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	O22850
Product Type	Recombinant Protein
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Purity	>85% (SDS-PAGE)
Sequence	TSKIKKFI LFLGVAVFY LYRYPSSPST VEQSSTSIYN ISVKDIEGKD VLSLKFTGKV LLIVNVASKC GLTHGNYKEM NILYAKYKTQ GFEILAFPCN QFGSQEPGSN MEIKETVCNI FKAEFPIFDK IEVNGKNTCP LYNFLKEQKG GLFGDAIKWN FAKFLVDRQG NVVDRYAPTT SPLEIEKDIV KLLASA
Source	Mammalian cell
Target Names	GPX3
Protein Names	Recommended name: Probable glutathione peroxidase 3, mitochondrial EC=1.11.1.9
Expression Region	13-206
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.