



# Recombinant Archaeoglobus fulgidus Geranylgeranylglyceryl phosphate synthase (AF\_0403)

<b>Product Code</b>	CSB-EP519905DOC-B
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	O29844
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Archaeoglobus fulgidus (strain ATCC 49558 / VC-16 / DSM 4304 / JCM 9628 / NBRC 100126)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MRWRKWRHIT KLDPDRTNTD EIIKAVADSG TDAVMISGTQ NVTYEKARTL IEKVSQYGLP IVVEPSDPSN VVYDVDYLFV PTVLNSADGD WITGKHAQWV RMHYENLQKF TEIIESEFIQ IEGYIVLNPD SAVARVTKAL CNIDKELAAS YALVGEKLFN LPIIYIEYSG TYGNPELVAE VKKVLDKARL FYGGGIDSRE KAREMLRYAD TIIVGNVIYE KGIDAFLETL P
<b>Source</b>	E.coli
<b>Target Names</b>	AF_0403
<b>Protein Names</b>	Recommended name: Geranylgeranylglyceryl phosphate synthase Short name= AfGGGPS Short name= GGGP synthase Short name= GGGPS EC= 2.5.1.41 Alternative name(s): (S)-3-O-geranylgeranylglyceryl phosphate synthase Phospho
<b>Expression Region</b>	1-231
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.