



Recombinant Staphylococcus aureus Glycyl-glycine endopeptidase lytM (lytM)

Product Code	CSB-BP520167FLF
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	O33599
Product Type	Recombinant Protein
Immunogen Species	Staphylococcus aureus (strain NCTC 8325)
Purity	>85% (SDS-PAGE)
Sequence	AETTN TQQAHTQMST QSQDVSYGTY YTIDSNNGDYH HTPDGNWNQA MFDNKEYSYT FVDAQGHTHY FYNCYPKNAN ANSGSQTYVN PATAGDNNDY TASQSQQHIN QYGYQSNVGP DASYYSHSNN NQAYNSHDGN GKVNYPNGTS NQNGGSASKA TASGHAKDAS WLTSRKQLQP YGQYHGGGAH YGVDYAMPEN SPVYSLTDGT VVQAGWSNYG GGNQVTIKEA NSNNYQWYMH NNRLTVSAGD KVKAGDQIAY SGSTGNSTAP HVHFQRMSGG IGNQYAVDPT SYLQSR
Source	Baculovirus
Target Names	lytM
Protein Names	Recommended name: Glycyl-glycine endopeptidase lytM EC= 3.4.24.75 Alternative name(s): Autolysin lytM
Expression Region	26-316
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.