



# Recombinant *Synechococcus elongatus* ATP-dependent Clp protease ATP-binding subunit ClpX (clpX)

|                          |   |
|--------------------------|---|
| <b>Product Code</b>      | CSB-YP520184FPY   |
| <b>Storage</b>           | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.   |
| <b>Uniprot No.</b>       | O34126  |
| <b>Product Type</b>      | Recombinant Protein   |
| <b>Immunogen Species</b> | <i>Synechococcus elongatus</i> (strain PCC 7942) ( <i>Anacystis nidulans</i> R2)  |
| <b>Purity</b>            | >85% (SDS-PAGE)   |
| <b>Sequence</b>          | MSRYDShLkC SFCGKSQEQV RkLIAGPGVY ICDECVDLCN EILDEELFDG<br>NANPAPPPGS PKPATAATTS HQGPRKPKRS LTLSQIPKPH DIKRYLDHHV<br>IGQNEAKKML SVAVYNHYKR LAYQQSAEEA GETIEIQKSN ILLIGPTGCG<br>KTLAQTLAK LLDVpFAVAD ATTLTEAGYV GEDVENILLR LLQVADLDIE<br>EAQRGIYID EIDKIARKSE NPSITRDVSG EGVQQALLKM LEGTIANVPP<br>QGGRKHPYQD CIQIDTRNIL FICGGAFVGL EKLVDrrLGK KsigfVhPEE<br>SKTKEQRAAA ILRHMEPEDL VQFGLIPEFI GRMPVTAVLD PLDETALTEI<br>LTEPRDALVK QYQKLMHMDS VELEFRQDAI EAIAREAFRR KTGARALRGI<br>VEEIMLEVMY ELPSRQDVTH CTITREmVEK RSSADVLLLp SSLPTPESA |
| <b>Source</b>            | Yeast   |
| <b>Target Names</b>      | clpX  |
| <b>Protein Names</b>     | Recommended name: ATP-dependent Clp protease ATP-binding subunit ClpX   |
| <b>Expression Region</b> | 1-449   |
| <b>Notes</b>             | Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.   |
| <b>Tag Info</b>          | Tag type will be determined during the manufacturing process.   |
| <b>Protein Length</b>    | full length protein   |
| <b>Reconstitution</b>    | We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.   |
| <b>Shelf Life</b>        | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.   |