



Recombinant *Schizosaccharomyces pombe* Uncharacterized but2-like protein C27D7.11c (SPAC27D7.11c)

Product Code	CSB-YP520369SXV
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	O42665
Product Type	Recombinant Protein
Immunogen Species	<i>Schizosaccharomyces pombe</i> (strain 972 / ATCC 24843) (Fission yeast)
Purity	>85% (SDS-PAGE)
Sequence	SPLRRRD DPLNGRQFGL LALHSGNEYV HLHGFYVGDS GSVYLDPADG TDDAATFTMS NGQLSVGTRY ASVSANGEIV FKSDTNSTSS GFSVGEAISS GYSLKYNTE SAVACPSLVN NQVYQVFFGV GNGNPSCVGI AVLAVLPQPI SSSSTYNSTT SSYHNSTSTP PPTITSTKAS TVTSTEATTV TTTTAVTVTA TETYTVTATN GGSTITSTGA STVTSTQPST VTSTQRKNTA TTTKTTTTYVG PSPSASSVVA YTTKCIVVPV ITTAASQSEA ATPSPSAAVY PLFPHGIRLI DSTKPEANSV NVYSPVVFQK QNNHTNTIFT FDVPQVSGSC ELNFHLDTSV FPITVEGVNG TGRFILFNLS SVANDSTVYS NRPNRIAEIG RFNCSSSGCD YATNVTCPDS YTAVSYEMMA LTDDSYLSFF EEADPLEGLT LRV
Source	Yeast
Target Names	SPAC27D7.11c
Protein Names	Recommended name: Uncharacterized but2-like protein C27D7.11c
Expression Region	24-463
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.