



# Recombinant *Trichophyton verrucosum* Probable aspartic-type endopeptidase CTSD (CTSD)

<b>Product Code</b>	CSB-MP520386TQC
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	D4D8U6
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	<i>Trichophyton verrucosum</i> (strain HKI 0517)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	FYPIKLPDF TKGLVSNHGS IDRRFFTFPG LYKHAHTGST TLNIRRGPSN YRRDNNYPAQ IASPTAPNT LGINNDGYDF SYFSEVKVGS EGQKMWMLID TGASGTWVFG SDCTSKACGR HNTFGKEDSK TIKVTDEKVG VTYGTGKVS VIVNDTMSFA GFELVTPFGS ASTASDDFLN YPMDGILGIG PQDPNAKTPT VVQLLMQQL LKSNVIGINL QRASEGATDG QITFGDIDKS KFSGELIYSN VVPDGYQWEI AMDDLIMDGK SLNLKGRGTGI IDTGTSFLIL PPADADLIHS MIPQANKGSG FYTLPCSTKV DIKLSIGGVE YTIQPDDYVG NETATKGTCN SLIVGRQILG PKQWLVDVDF LKNVYSVDFD DKNRVGLAAR KYAGTKNPPS STPSPGMFLL HAILCPKTIS VLMLHIDPTS NKAPSGGSPG LPAESGSDST TNGEATNGAT SSPNS
<b>Source</b>	Mammalian cell
<b>Target Names</b>	CTSD
<b>Protein Names</b>	Recommended name: Probable aspartic-type endopeptidase CTSD EC=3.4.23.-
<b>Expression Region</b>	22-485
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	Full Length of Mature Protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.