



Recombinant *Saccharomyces cerevisiae* Protein HRI1 (HRI1)

Product Code	CSB-EP520685SVO
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	E7KS00
Product Type	Recombinant Protein
Immunogen Species	<i>Saccharomyces cerevisiae</i> (strain Lalvin QA23) (Baker's yeast)
Purity	≥85% (SDS-PAGE)
Sequence	MPALLKRLLF QVGPHPNERT FTLSSVSTDG HYISLRPFVK PSGDELSFPF EWAFAGTNET VKVNDQNGV VTQDFNFWLD TNVYLNVPNT HRGEVNTTWK NWDSGCVEET GAVYPPGADK ESVSFRELWQ PVDPSREDLV IVSPNBEKFS SNARSIVLKV TDEAYDGLVI VIGRWIQGFL SQKNNNTIEG LNFIRLLEKD SGKSEFLLSY GKEVNIKIPQS YENLKKGSTV TSNGLNWEVI EYHA
Source	E.coli
Target Names	HRI1
Protein Names	Recommended name: Protein HRI1 Alternative name(s): HRR25-interacting protein 1
Expression Region	1-244
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.