



# Recombinant Human RNA 3'-terminal phosphate cyclase (RTCD1)

<b>Product Code</b>	CSB-EP520755HU-B
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	O00442
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Purity</b>	≥85% (SDS-PAGE)
<b>Sequence</b>	MAGPRVEVDG SIMEGGGQIL RVSTALSCLL GLPLRVQKIR AGRSTPGLRP QHLSGLEMIR DLCDGQLEGA EIGSTEITFT PEKIKGGIHT ADTKTAGSVC LLMQVSMPCV LFAASPSELH LKGGTNAEMA PQIDYTMVF KPIVEKFGFI FNCDIKTRGY YPKGGGEVIV RMSPVKQLNP INLTERGCVT KIYGRAFVAG VLPFKVAKDM AAAAVRCIRK EIRDLVNIQ PVQEPKDQAF GNGNGIIIA ETSTGCLFAG SSLGKRGVNA DKVGIEAAEM LLANLRHGGT VDEYLQDQLI VFMALANGVS RIKTGPVTLH TQTAIHFAEQ IAKAKFIVKK SEDEEDAAKD TYIIECQGIG MTNPNL
<b>Source</b>	E.coli
<b>Target Names</b>	RTCA
<b>Protein Names</b>	Recommended name: RNA 3'-terminal phosphate cyclase Short name= RNA cyclase Short name= RNA-3'-phosphate cyclase EC= 6.5.1.4 Alternative name(s): RNA terminal phosphate cyclase domain-containing protein 1 Short name= RT
<b>Expression Region</b>	1-366
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.