



Recombinant Shigella flexneri tRNA(fMet)-specific endonuclease VapC (vapC)

Product Code	CSB-MP520879SZB
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	O06662
Storage Buffer	Lyophilized from Tris/PBS-based buffer, 6% Trehalose, pH 8.0
Product Type	Recombinant Proteins
Immunogen Species	Shigella flexneri
Purity	>85% (SDS-PAGE)
Sequence	MLKFMLDTNI CIFTIKNKPA SVRERFNLNQ GKMCISSVTL MELIYGAEKS QMPERNLAVI EGFVSRIDVL DYDAAAATHT GQIRAEARQ GRPVGPFQDM IAGHARSRL IIVTNTREF ERVGGLRTED WS
Source	Mammalian cell
Target Names	vapC
Protein Names	Recommended name: tRNA(fMet)-specific endonuclease VapC EC= 3.1.-.- Alternative name(s): RNase VapC Toxin VapC
Expression Region	1-132
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.