



Recombinant Enterobacter aerogenes Uridine phosphorylase (udp)

Product Code	CSB-EP520955EKN-B
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	O08444
Product Type	Recombinant Protein
Immunogen Species	Klebsiella aerogenes (Enterobacter aerogenes)
Purity	>85% (SDS-PAGE)
Sequence	MSKSDVFHLG LTKNDLQGAT LAIVPGDPER VEKIAALMDK PVKLASHREF TSWRAELDGK PVIVCSTGIG GPSTSIAVEE LAHVGVRTFL RIGTTGAIQP HINVGDLVLT TGSVRLDGAS LHFAPMEFPA VADFACTTAL VEAAKSIGAT THIGVTASSD TFYPGQERYD TFSGRVVSFR KGSMEEWQAM GVMNYEMESA TLLTMCASQG LRAGMVAGVI VNRTQQEIPN AETMKQTESH AVKIVVEAAR RLL
Source	E.coli
Target Names	udp
Protein Names	Recommended name: Uridine phosphorylase Short name= UPase Short name= UrdPase EC= 2.4.2.3
Expression Region	1-253
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.