



Recombinant *Schizosaccharomyces pombe* Serine/threonine-protein kinase haspin homolog C23C4.03 (SPAC23C4.03)

Product Code	CSB-BP521064SXV
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	O13924
Product Type	Recombinant Protein
Immunogen Species	<i>Schizosaccharomyces pombe</i> (strain 972 / ATCC 24843) (Fission yeast)
Purity	>85% (SDS-PAGE)
Sequence	MESRSKVKTY GKNRRYINKD IEIWASLDER PCALKPRNIE NFNNQERSDL EHIHSPKKD SLLSWNILLK KGSYKENELL AKRNQNLVPT VIIPASPRDN ASKSVVSKKE VVNLSSSVL SGKPANNSKL DPLHRLLQIV AQEDALPFSQ FVKSQTFEIQ KIGEASYSEV YQASNADDVP VVWKVIPFGE DGQAQYADVL NEVQISQWIK VDGFANLHQV VVVKGTYP SL LLEEWDRYLM QNGSENDERPD SYSSTQLYCV LCLDHSGTDL EHFELRSWRE CWSV FYETLK ILSLVETRYE FEHRDLHWGN ILIRKADRSE EEVSFLLNEI SLDDIESVDF PGSQDKADDF DNILQVTLID FTLARASY SQ GIISYNEFND PDLFNGVDDY QFDIYRLMSR VTKGRWAQFF PITNVLWLHY LIHQLLHKKN LSSPLTETET LMRSRLKQIF RLIDPVKTMQ FQQAEDSIRS KSTVTSATSL LNWVRQKY
Source	Baculovirus
Target Names	hrk1
Protein Names	Recommended name: Serine/threonine-protein kinase haspin homolog C23C4.03 EC= 2.7.11.1
Expression Region	1-488
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.



Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.