



Recombinant *Schizosaccharomyces pombe* Serine--tRNA ligase, cytoplasmic (SPAC29A4.15)

Product Code	CSB-MP521079SXV
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	O14018
Product Type	Recombinant Protein
Immunogen Species	<i>Schizosaccharomyces pombe</i> (strain 972 / ATCC 24843) (Fission yeast)
Purity	≥85% (SDS-PAGE)
Sequence	MLDINLFQVE KGGNPEIIRE SQRKRGADVG VVDKVIEMYK EWVSLRFELD NTNKSINRVQ KEIGLKMKAK EDASELLEEK NSLTERKKNL IEQETAKNKE MLNVVSSIGN IVHDSVPVSM DEDNNEIIRK WAPEGVTVEK KNCLSHHEVL TRLDGYDPER GVKVSGHRGY FLRQYGVFFN LALIYGLDF LEKRGYIALQ APTMLNKDVM AKTAQLEQFD EELYKVIDGD EERYLIATSE QPISAYHSGE WFEKPSEQLP LKYAGYSTCY RREAGSHGRD AWGIFRVHAF EKIEQFVLTD PEKSWEAFTE MINHAEDFYK SLELPYRIVA IVSGALNNA AKKYDLEAWF PFQGEYKELV SCSNCTDYQS RNLEIRCGVK KMGDREKKYV HCLNSTLCAT ERALCCILEN YQTPDGVNVP KVLQPYMG GK TFLPFTKELP KNSTSKKGKN
Source	Mammalian cell
Target Names	SPAC29A4.15
Protein Names	Recommended name: Serine--tRNA ligase, cytoplasmic EC= 6.1.1.11 Alternative name(s): Seryl-tRNA synthetase Short name= SerRS Seryl-tRNA(Ser/Sec) synthetase
Expression Region	1-450
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.