



Recombinant *Schizosaccharomyces pombe* Putative methyltransferase C2C4.06c (SPAC2C4.06c)

Product Code	CSB-BP521083SXV
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	O14039
Product Type	Recombinant Protein
Immunogen Species	<i>Schizosaccharomyces pombe</i> (strain 972 / ATCC 24843) (Fission yeast)
Purity	>85% (SDS-PAGE)
Sequence	MDFYNHAANI LSDLSKKKGS IKQLAFNSKK HDPKRTYALV CETLKYKPVL DEIARSELL VLEKCLKENL ARVLVHLLM SKRGLSISNG PIKECILRHK TRLNAEFVKL KVKKGVKSHE ELALKNPVSL PRWLRINTIK STKDEVLQGL GLDKVSSIEE LGPDKFYIDD CVENLIAIDP SFPIVENSLY KEGKVIIQDK ASCFPAAVLA GLTGHVGDII DGCAAPGNKT THLAACFPKS HIFAFERDAK RVQTLRKMVG ISGANNVTIE HQDFTLTDPK SDLYRNVTHI LLDPSCSGSG IVSRQDYLLG NEQDVTEDETE RLENLCSFQS TILKHALQFP NCRHVTYSTC SVHRENEQV VCEVLSQEPD WKCNSLTKTL PNWKTRGIPE YCAQPSMAEG MIRCKPGAGG TIGFFVANLY HPQREQETFK MYKNDDDTKK RKRKKKKKKEV KKKARIQGEE
Source	Baculovirus
Target Names	rcm1
Protein Names	Recommended name: Putative methyltransferase C2C4.06c EC= 2.1.1.-
Expression Region	1-460
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.