



Recombinant Schizosaccharomyces pombe Ribosomal RNA-processing protein 14-N (rrp14n)

Product Code	CSB-EP521123SXV
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	O14279
Product Type	Recombinant Protein
Immunogen Species	Schizosaccharomyces pombe (strain 972 / ATCC 24843) (Fission yeast)
Purity	>85% (SDS-PAGE)
Sequence	MTEVEERINA WNDGFEELLS YIPAKLYYRE HTANQWKQKK STLEEKKERR KMKFSLDGLV NDADDDSQKS SKWEQSSTMS DTDVSDRHA ESMSQIRGKL ASKIQLREK RKAGDLNQKR QNKRPVENEK DSQKGGKSK VQKKKHKASP RAGF
Source	E.coli
Target Names	rrp14n
Protein Names	Recommended name: Ribosomal RNA-processing protein 14-N Alternative name(s): Ribosome biogenesis protein rrp14-N
Expression Region	1-154
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.