



Recombinant Uncharacterized protein K07E8.7 (K07E8.7)

Product Code	CSB-MP521167CXY
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	O16686
Product Type	Recombinant Protein
Immunogen Species	Caenorhabditis elegans
Purity	>85% (SDS-PAGE)
Sequence	MSDTTDPEN QKSPKPSGKA DKRKIEEKPE NSSLKRKKFE DPNKKVDPLE ELPMKVPFKI VDGVRHLAPY WACYRTRTKG RWIGRKMVEV FSGEFLSTNR NYAKIACKMG RIYVNGEQMT DVDYVMRNGD RVEHWahrHE HPIRDLPIRV ISETDDLFFV EKPPSLPVHT CGQYAIHTVL GQLRVNEGRT GLRVLHRLDR ATSGVLLFAK NYETDLEFKT TLKQGEWSKE YICKVDGVFP DEEQVCEQPI GPLVISMGIQ CVRPDGDGKDAK SRFRKLWSDG TQSVVQVHIE TGRTHQIRVH SQFLGHPIAG DQIYNSAVWG PTKGKNADYQ KSFDELCEDV RNTHKCNWH EKPNPEFEQR MEHLAADTTP ITPEAPSLTL EQRPEFDEIC QKCNVESKKV PENHFQLYLH CLKYETKKWS FKTEMPDWAV QK
Source	Mammalian cell
Target Names	K07E8.7
Protein Names	Recommended name: Uncharacterized protein K07E8.7
Expression Region	1-432
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.