



Recombinant Putative H/ACA ribonucleoprotein complex subunit 4 (K01G5.5)

Product Code	CSB-BP521186CXY
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	O17919
Product Type	Recombinant Protein
Immunogen Species	Caenorhabditis elegans
Purity	≥85% (SDS-PAGE)
Sequence	MGKKDKRSLK EGDDLAEAAQK KGSFQLPSSN ETAKLDASQW PLLLNKYDKL NVRTNHYTPH VEGVSPVKRD IKNYISSGFF NLDKPSNPSS HEVVSWKRI LRCEKTGHSG TLDPKVSGCL IVCIDRTTRL AKSQQGAGKE YICIFKLHEE VEDDRKVKQA LEKLTGALFQ RPPLISAVKR QLRIRTVYEN KFIEYDPAQQ MGIFNCICES GTYVRTICVH LGLILGCGGQ MQELRRNRSG ICDENENMVT MHDVLDQAQYL LDTQKDESYM RHIVRPLEAL LTQHKRVVVK DSCINAICYG AKILIPGILR YDDDIEVGKE IVIMSTKGEA ICIAIAQMNT STIASVDHGV VAKSKRVIME RDVYGRKWGL GPVASKKKQM VKDGLLDKFG KPNDTTPKSW AKEYVQTSTK KEVKKEETPD EEEEEAPKKK SKSKKKQESS DSDSD
Source	Baculovirus
Target Names	K01G5.5
Protein Names	Recommended name: Putative H/ACA ribonucleoprotein complex subunit 4 EC=5.4.99.-
Expression Region	1-445
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.