



Recombinant Uncharacterized protein ZK849.5 (ZK849.5)

Product Code	CSB-BP521191CXY
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	O18304
Product Type	Recombinant Protein
Immunogen Species	Caenorhabditis elegans
Purity	≥85% (SDS-PAGE)
Sequence	MTVPYYRDVP TGNALRNFFK IMLKRKGNLF TAIFKELCLF LGLYFLFMVI YRLVLPKSGG QREIRKIVEN LLSHQEFTIP LEFLLGFFVS SVVDRWRKSL DSLAYIESCA HAVVIGFPPN SNGSDKNLLI RRTIIRYLVV SQILLYREIS LKVRRRFKKL TILGKAKLLN QNEIDKLNKL ECKHYDIYFL PISWALSIE DKIDKENLAN EFTILWGQIK EWQTKLSLLR NCDYIPIPLA YPQAVFLAVR CYFAVCVFTR QHLDYDYSKM HTWITFFPVL TTFQYIFMMG WMKVAEILLN PMGEDEDDFE LNFIIDNNLK NGLDIVSGLC GNHRKLAEHE IENDCRPYYQ TNEQDRKKNR APPELKNVE FKSFTMEKAS KDSPMVVKSE SVCCLRRRFN SSAKKHHDSK V
Source	Baculovirus
Target Names	best-26
Protein Names	Recommended name: Uncharacterized protein ZK849.5
Expression Region	1-411
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.