



# Recombinant Arabidopsis thaliana NADPH-dependent thioredoxin reductase 3 (NTRC)

<b>Product Code</b>	CSB-EP521234DOA-B
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	O22229
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Arabidopsis thaliana (Mouse-ear cress)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	ATA NSPSSSSSGG EIIENVVIIG SGPAGYTAAI YAARANLKPV VFEGYQMGGV PGGQLMTTTE VENFPGFDPG ITGPDLMEKM RKQAERWGAE LYPEDVESLS VTTAPFTVQT SERKVKCHSI IYATGATARR LRLPREEEFW SRGISACAIC DGASPLFKGQ VLAVVGGGDT ATEEALYLT YARHVHLLVR RDQLRASKAM QDRVINNPNI TVHYNTETVD VLSNTKGQMS GILLRRLDTG EETELEAKGL FYGIGHSPNS QLLEGQVELD SSGYVLVREG TSNTSVEGVF AAGDVQDHEW RQAVTAAGSG CIAALSAERY LTSNNLLVEF HQPQTEEAKK EFTQRDVQEK FDITLTKHKG QYALRKLYHE SPRVILVLYT SPTCGPCRTL KPILNKVVDE YNHVHFVEI DIEEDQEIAE AAGIMGTPCV QFFKNKEMLR TISGVKMKKE YREFIEANK
<b>Source</b>	E.coli
<b>Target Names</b>	NTRC
<b>Protein Names</b>	Recommended name: NADPH-dependent thioredoxin reductase 3 Short name= NTR3 EC= 1.8.1.9 Alternative name(s): NADPH-dependent thioredoxin reductase C Short name= ANTR-C Short name= AtNTRC
<b>Expression Region</b>	68-529
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	Full Length of Mature Protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.