



Recombinant Arabidopsis thaliana Peroxidase 19 (PER19)

Product Code	CSB-BP521262DOA
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	O22959
Product Type	Recombinant Protein
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Purity	>85% (SDS-PAGE)
Sequence	KPPAPRPHR ELSADYYSKK CPQLETLVGS VTSQRFKEVP ISAPATIRLF FHDCFVEGCD GSILIEKKG SKKLAEREAY ENKELREEGF DSIKAKALV ESHCPSLVSC SDILAI AARD FIHLAGGPYY QVKKGRWDGK RSTAKNVPPN IPRSNSTVDQ LIKLFASKGL TVEELVVLGS SHTIGFAHCK NFLGRLYDYK GTKRPDPSLD QRLLKELRMS CPFSGGSSGV VLPLDATTFF VFDNGYFTGL GTNMGLLGSD QALFLDPRTK PIALEMARDK QKFLKAFGDA MDKMGSIGVK RGKRHGEIRT DCRVFL
Source	Baculovirus
Target Names	PER19
Protein Names	Recommended name: Peroxidase 19 Short name= Atperox P19 EC= 1.11.1.7 Alternative name(s): ATP51
Expression Region	32-346
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.