



Recombinant Arabidopsis thaliana Putative F-box/LRR-repeat protein At4g13960 (At4g13960)

Product Code	CSB-YP521272DOA
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	O23257
Product Type	Recombinant Protein
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Purity	≥85% (SDS-PAGE)
Sequence	MDHVSSLPDE VLYHILSFLT TKEAALTSIL SKRWRNLFTF VPNLIDIDDSV FLHPQEGKED RYEIQKSFMK FVDRVLALQG NSPIKKLSLK LRTGFDSHRV DGWISNALAR GVTELDLLII LNLVTLKLNS VRVDWLAAGD IFLPMLKTLV LHSVRLCVDN VNWGGKDVNV SNASLKILTV NYNICLGTFVS VDTPSLVYFC FTEYVAIDYP VVIMQNLFEA RISLLVTQDQ IERARAPNID WVEDDVVLRL VNMEKLIKGI RNVQYMYLSP NTLEVLSLCC ESMPVFNNLK SLTLKSHKSR GWQAMPVLLR NCPHLETLLV EGILHHVTDK CGYACDCVSR EEKGCSLTSC PVKVLEIKGF RGTIKEMHMI KHFNYYFPCL METNIYIEES GPPKLRVPQV CEFIAEMMEH YGKSSSFNVQ LLVGGYLYKK WTAQ
Source	Yeast
Target Names	At4g13960
Protein Names	Recommended name: Putative F-box/LRR-repeat protein At4g13960
Expression Region	1-434
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.