



Recombinant Arabidopsis thaliana Proteasome subunit beta type-7-A (PBB1)

Product Code	CSB-EP521284DOA
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	O23710
Product Type	Recombinant Protein
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Purity	>85% (SDS-PAGE)
Sequence	TGT TIVGLIFKDG VILGADTRAT EGPIVADKNC EKIHYPAPNI YCCGAGTAAD TEAVTDMVSS QLRLHRYQTG RDSRVITALT LLKKHLFSYQ GHVSAALVLG GVDITGPHLH TIYPHGSTD T LPFATMGSGS LAAMSVFEAK YKEGLTRDEG IKLVAESICS GIFNDLGSGS NVDICVITKG NKEYLRNYME PNPRTYVSSK GYSFTKKTEV LLTKITPLLE RVEITEVGEA MEE
Source	E.coli
Target Names	PBB1
Protein Names	Recommended name: Proteasome subunit beta type-7-A EC= 3.4.25.1 Alternative name(s): 20S proteasome beta subunit B-1 Proteasome component FA Proteasome component FB Proteasome subunit beta type-2
Expression Region	38-273
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.