



Recombinant Arabidopsis thaliana Hydroxyacylglutathione hydrolase cytoplasmic (GLX2-2)

Product Code	CSB-MP521308DOA
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	O24496
Product Type	Recombinant Protein
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Purity	>85% (SDS-PAGE)
Sequence	MKIFHVPCLQ DNYSYLIIDE STGDAAVVDP VDPEKVIASA EKHQAKIKFV LTTHHHWDHA GGNEKIKQLV PDIKVYGGSL DKVKGCTDAV DNGDKLTLGQ DINILALHTP CHTKGHISYY VNGKEGENPA VFTGDTLFVA GCGKFFEGTA EQMYQSLCVT LAALPKPTQV YCGHEYTVKN LEFALTVEPN NGKIQQKLAW ARQQRQADLP TIPSTLEEEL ETNPFMRVVK PEIQEKLKCK SPIDTMREVR NKKDQWRG
Source	Mammalian cell
Target Names	GLX2-2
Protein Names	Recommended name: Hydroxyacylglutathione hydrolase cytoplasmic EC= 3.1.2.6 Alternative name(s): Glyoxalase II Short name= Glx II
Expression Region	1-258
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.