



Recombinant Zea mays CLAVATA3/ESR (CLE)-related protein 2-B

Product Code	CSB-EP521314ZAX
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	O24618
Product Type	Recombinant Protein
Immunogen Species	Zea mays (Maize)
Purity	>85% (SDS-PAGE)
Sequence	NVWQ TDEDAFYSTN KLGVNGNMEM AQQQGGFIGH RPRLASFNRA SKQLDREKRP VPSGPDPIHH SIPSHAPQHP PSYGKAPYED DKSIASPLS NLIGPPPFLD RY
Source	E.coli
Target Names	ESR2B
Protein Names	Recommended name: CLAVATA3/ESR (CLE)-related protein 2-B Alternative name(s): Embryo surrounding region protein 2-B Cleaved into the following chain: 1. ESR2Bp
Expression Region	27-132
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.