



# Recombinant Bacillus subtilis Glutaminase 1 (glsA1)

<b>Product Code</b>	CSB-MP521711BRJ
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	O31465
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Bacillus subtilis (strain 168)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MKELIKEHQK DINPALQLHD WVEYYRPF AA NGQSANYIPA LGKVNDSQLG ICVLEPDGTM IHAGDWNVSF TMQSISKVIS FIAACMSRGI PYVLDRVDVE PTGDAFNSII RLEINKPGKP FNPMINAGAL TIASILPGES AYEKLEFLYS VMETLIGKRP RIHEEVFRSE WETAHRNRAL AYYLKETNFL EAEVEETLEV YLKQCAMEST TEDIALIGLI LAHDGYHPIR HEQVIPKDVA KLAKALMLTC GMYNASGKYA AFVGVPKSG VSGGIMALVP PSARREQPFQ SGCIGIGIYGP AIDEYGNSLT GGMLLKHMAQ EWELSIF
<b>Source</b>	Mammalian cell
<b>Target Names</b>	glsA1
<b>Protein Names</b>	Recommended name: Glutaminase 1 EC= 3.5.1.2
<b>Expression Region</b>	1-327
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.