



Recombinant Homospermidine synthase (hss)

Product Code	CSB-BP521859RLX
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	O32323
Product Type	Recombinant Protein
Immunogen Species	Blastochloris viridis (Rhodospseudomonas viridis)
Purity	>85% (SDS-PAGE)
Sequence	MTDWPVYHRI DGPIVMIGFG SIGRGTLP LI ERHFAFDRSK LVVIDPSDEA RKLAEARGVR FIIQAVTRDN YRELLVPLLT AGPGQGFCVN LSVDTSSLDI MELARENGAL YIDTVVEPWL GFYFDPDLKP EARSNYALRE TVLAARRNKP GGTTAVSCCG ANPGMVSWFV KQALVNLAAD LGVTGEEPTT REEWARLAMD LGVKGIIHAE RDTQRASFPK PFDVFVNTWS VEGFVSEGLQ PAELGWGTFE RWMPDNARGH DSGCGAGIYL LQPGANTRVR SWTPTAMAQY GFLVTHNESI SIADFLTVRD AAGQAVYRPT CHYAYHPCND AVLSLHEMFG SGKRQSDWRI LDETEIVDGI DELGVLLYGH GKNAYWYGSQ LSIEETRIA PDQATGLQV SSAVLGAMVW ALENPNAGIV EADDLDFRRC LEVQTPYLG P VVG VYTDWTP LAGRPGLFPE DIDTSDPWQF RNVLV RD
Source	Baculovirus
Target Names	hss
Protein Names	Recommended name: Homospermidine synthase Short name= HSS EC= 2.5.1.44
Expression Region	1-477
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.