



Recombinant *Schizosaccharomyces pombe* Putative methyltransferase-like protein C27D7.08c (SPAC27D7.08c)

Product Code	CSB-EP522112SXV-B
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	O42662
Product Type	Recombinant Protein
Immunogen Species	<i>Schizosaccharomyces pombe</i> (strain 972 / ATCC 24843) (Fission yeast)
Purity	>85% (SDS-PAGE)
Sequence	MPIYILERS VKNGRIDFWN EDAIRTLGKA ILDRDYSLRV EFPENRLCPM VPNRATYIRY IHDLLSSTSG QKDKKRIIGL DIGTGASCIY PLLGCRMYSY DFVGTEIDKF SFETAKSNIL QNNMESQIKI VLRSKQDCLL PDTEGMEEFT FVMCNPPFYE HEEDFINFKQ NPPSGVCTGV YHEMVTEGGE VGFANKILTE SKKRKGIQWY TCMFGKKSSV PAVVDKLRQ NISNYGIYEL ALGKTKRWII CWSFQAMRPH NELIRPSSTS LSKYFPHKVL QNWTLDPCLC AQIDDILQKF LDDNKIPWSK KGSVLEISTK SITWSRKARR ISKSQTSVSS LEGQMKCELN VIDNQLQCKW IEGYDYNVYE SFCSALARAL RDNKK
Source	<i>E.coli</i>
Target Names	SPAC27D7.08c
Protein Names	Recommended name: Putative methyltransferase-like protein C27D7.08c EC=2.1.1.-
Expression Region	1-385
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.