



Recombinant Agaricus bisporus Polyphenol oxidase 2 (PPO2)

Product Code	CSB-MP522115AYJ
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	O42713
Product Type	Recombinant Protein
Immunogen Species	Agaricus bisporus (White button mushroom)
Purity	≥85% (SDS-PAGE)
Sequence	MSLIATVGPT GGVKNRLNIV DFKNEKFFT LYVRSLELLQ AKEQHDYSSF FQLAGIHGLP FTEWAKERPS MNLYKAGYCT HGQVLFPTWH RTYLSVLEQI LQGAAIEVAK KFTSNQTDWV QAAQDLRQPY WDWGFELMPP DEVIKNEEVN ITNYDGKKIS VKNPILRYHF HPIDPSFKPY GDFATWRRTV RNPDRNRRED IPGLIKMRL EEGQIREKTY NMLKFNDawe RFSNHGISDD QHANSLESVH DDIHMVVGYG KIEGHMDHPF FAAFDPIFWL HHTNVDRLLS LWKAINPDVW VTSGRNRDGT MGIAPNAQIN SETPLEPFYQ SGDKVWTSAS LADTARLGYS YPDFDKLVGG TKELIRDAID DLIDERYG
Source	Mammalian cell
Target Names	PPO2
Protein Names	Recommended name: Polyphenol oxidase 2 Short name= PPO2 Short name= Phenolase 2 EC= 1.14.18.1 Alternative name(s): Cresolase 2 Tyrosinase 2
Expression Region	1-378
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.