



Recombinant *Saccharomyces cerevisiae* Suppressor of hydroxyurea sensitivity protein 2 (SHU2)

Product Code	CSB-YP522248SVR
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	E7QCS8
Product Type	Recombinant Protein
Immunogen Species	<i>Saccharomyces cerevisiae</i> (strain Zymaflore VL3) (Baker's yeast)
Purity	≥85% (SDS-PAGE)
Sequence	MSKDVIEYSK LFAKLVNTND DTKLDDTIAS FLYYMFREL FIRAISSL SDMFIYILDR VHNKEGNEHT SLIDVLVDEF YKGSSNSLLE YRLIVKDTND GAPPILVDIA HWFCSCEEFC KYFHEALEKT DEKEELHDVL INEVDDHLQF SDDRFAQLDP HSLSKQWYFK FDKICCSHLL AFSILLRSSI NVLKFFTVNS NKVFVIAIDN IDEWLNHIN IVE
Source	Yeast
Target Names	SHU2
Protein Names	Recommended name: Suppressor of hydroxyurea sensitivity protein 2
Expression Region	1-223
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.