



# Recombinant Arabidopsis thaliana Polyadenylate-binding protein RBP47A (RBP47A)

<b>Product Code</b>	CSB-BP522272DOA
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	F4I3B3
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Arabidopsis thaliana (Mouse-ear cress)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	<p>MQTPNNGST DSVLPPTSAG TTPPPPLQQS TPPPQQQQQQ            QWQQQQQWMA AMQQYPAAAM AMMQQQQMMM YPHPQYAPYN            QAAYQQHPQF QYAAAYQQQQQ QHHQSQQQPR GGSGGDDVKT            LWVGDLLHWM DETYLHTCFS HTNEVSSVKV IRNKQTCQSE GYGFVEFLSR            SAAEEALQSF SGVTMPNAEQ PFRLNWASFS TGEKRASENG PDLISFVGDL            APDVSDAVLL ETFAGRYPSV KGAKVVIDSN TGRSKGYGFV RFGDENERSR            AMTEMNGAFC SSRQMRVIA TPKRAAAYGQ QNGSQALTLA            GGHGGNGSMS DGESNNSTIF VGGLDADVTE EDLMQPFSDF GEVVSVKIPV            GKGCQFVQFA NRQSAEEAIG NLNGTVIGKN TVRLSWGRSP            NKQWRSDSGN QWNGGYSRGQ GYNNGYANQD SNMYATAAAA VPGAS</p>
<b>Source</b>	Baculovirus
<b>Target Names</b>	RBP47A
<b>Protein Names</b>	Recommended name: Polyadenylate-binding protein RBP47A Short name= Poly(A)-binding protein RBP47A Alternative name(s): RNA-binding protein 47A Short name= AtRBP47A
<b>Expression Region</b>	1-445
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.