



Recombinant Probable Hsp90 co-chaperone cdc37 (cdc-37)

Product Code	CSB-BP522327CXY
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	O02108
Product Type	Recombinant Protein
Immunogen Species	Caenorhabditis elegans
Purity	≥85% (SDS-PAGE)
Sequence	MPIDYSKWKD IEVSDDEDDT HPNIDTPSLF RWRHQARLER MAEKKMEQEK IDKEKGTTSK KMEELEKKLA AADVTDKSDI QKQIDEVKAQ EEAWRKKEAE LEEKERLEPW NVDTIGHEAF STSRINKITE KKPQAPKTDE EDTHAMSTFF ETHESLLEKM AVLKNGAKST ELFLAEHPHM ASEYTANWLT IEALNAAIDF NEEKMKTMAE QCIIIQYLLE LSKSLNAVAT NTTVQKQFFK KFEAAEPVYM KH YQDEVKAF EDRLRTRAQT KRDAAMEEAE AEEKAERMKS APGGIDPQEV FEQLPEEMRK CFEAHDIEAL KGVAQKMDEE VF KYHFDRCI ASGLWVPGKA DDDDDDDEEAA PAAAAEPTTSS
Source	Baculovirus
Target Names	cdc-37
Protein Names	Recommended name: Probable Hsp90 co-chaperone cdc37 Alternative name(s): Cell division cycle-related protein 37 Hsp90 chaperone protein kinase-targeting subunit
Expression Region	1-370
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.