



Recombinant AT-rich interactive domain-containing protein cfi-1 (cfi-1)

Product Code	CSB-EP522330CXY-B
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	O02326
Product Type	Recombinant Protein
Immunogen Species	Caenorhabditis elegans
Purity	≥85% (SDS-PAGE)
Sequence	MSVRIDEPQL FVSMSKEPTQ ETVNVGGHHD DSSSNCDERV DDQTEEQKSP PASPDLTANL NVFDLESRQK VVQRLLNSQL NLSNLRAPLN LPPIFQALQG PFSIQQLLG LASGLTAISP GLDDYDEENT NQGEPEDLTL GGFRKETSVK SEEPSSEGIN ASGPAWSYEE QFKQLYELSD DVKRKEWLDD WLNFMHRIGK PVTRIPIMAK QVLDLYELYR LVVQHGGGLVE IINKKLWREI TKGLNLPSSI TAAFTLRTQ YQKYLYDYEC EKEKLSNQSD LQQAIDGNRR EAPGRRTAPS FPLPFQLPHA ASAAATMLNN QLNLGLMRND LLDDENTLSL QASGLFGTSY GAEQMAILEA HQRNLERAQR AVQQQVARQS LGLTACSNGN GGNIHNSGRE STSSNDSDIP AKRPKLENDV KTINGASSMRI STKHSDNSKT SMSVSMEING ITYQGVLFAL DETVSES
Source	E.coli
Target Names	cfi-1
Protein Names	Recommended name: AT-rich interactive domain-containing protein cfi-1 Alternative name(s): ARID domain-containing protein CFI-1
Expression Region	1-467
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.