



Recombinant Arabidopsis thaliana Putative F-box/kelch-repeat protein At1g62270 (At1g62270)

Product Code	CSB-MP522359DOA
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	O04591
Product Type	Recombinant Protein
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Purity	>85% (SDS-PAGE)
Sequence	MTTLKTSSFS SLPWDLVEDI LARVPATSLK RLRSTCKQWN FLFNDQIFTK MHFDKAEKQF LVLILRLYTV CSMSLDLRGL HDNIDPSIEV KGELSLIDPH CSSRKTFVSK VFHCNGLLLC TTMTGLVWWN PCTDQTRWIK TEVPHNRNDK YALGYGNYKS CYNKIMKFL DLESFDLEIY EVNSNSWRVL GTVTPDFTIP LDAEGVSLRG NSYWIASHKR EEEEEEEEE NEYFINDFLI SFDFITTERFG PRVSLPFKCE SSWDTISLSC VREERLSLFF QDDGTLKMEI WMTNNITETK TTTMSWSPFL KIDLTYGHR FGNEVSFLVD EENKVIVCCD EEEDDINDTV YIIGENEYWR KEDIVQRSYR PRMFSYVPSL VQI
Source	Mammalian cell
Target Names	At1g62270
Protein Names	Recommended name: Putative F-box/kelch-repeat protein At1g62270
Expression Region	1-383
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.